

Product Datasheet

NUBP1 antibody (orb537187)

Description Rabbit polyclonal antibody to NUBP1

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Recombinant fusion protein of human NUBP1(NP_002475.2).

Preservatives PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 38kDa


Uniprot ID **P53384**

Dilution Range WB: 1:500-1:2000


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of mouse li...



Immunohistochemical staining of human co...



Immunohistochemical staining of human th...