



Product Datasheet

SERPINA10 antibody (orb537194)


Description	Rabbit polyclonal antibody to SERPINA10
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human SERPINA10 (NP_001094077.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	51kDa
Uniprot ID	Q9UK55
Entrez	51156
Dilution Range	WB: 1:500-1:2000, IF: 1:20-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of PC-12 cell...



Immunofluorescence analysis of A549 cell...



Western blot analysis of extracts of Hep...