


Product Datasheet

TNFRSF12A antibody (orb537218)

Description	Rabbit polyclonal antibody to TNFRSF12A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-129 of human TWEAKR/Fn14/CD266 (NP_057723.1).
Preservatives	PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14kDa
Uniprot ID	Q9NP84
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of lysates from HT...