

Product Datasheet

TNFRSF6B antibody (orb537245)

Description Rabbit polyclonal antibody to TNFRSF6B

Reactivity Human, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Recombinant fusion protein of human TNFRSF6B(NP_003814.1).

Preservatives PBS with 0.01% thiomersal,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 33KDa


Uniprot ID **O95407**

Dilution Range WB: 1:500-1:2000, IHC-P: 1:50-1:200


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of rat kidn...



Immunohistochemical staining of human ki...



Immunohistochemical staining of human st...