


Product Datasheet

VEGFR2 (Phospho-Y996) antibody (orb537270)

Description	Rabbit polyclonal antibody to VEGFR2 (Phospho-Y996)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y996 of human VEGF Receptor 2(NP_002244.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	200kDa
Uniprot ID	P35968
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HUV...