

Product Datasheet

PLCG1 (Phospho-S1248) antibody (orb537300)

Description Rabbit polyclonal antibody to PLCG1 (Phospho-S1248)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen A synthetic phosphorylated peptide around S1248 of human PLC gamma 1 (PLCG1) (NP_002651.2).

Preservatives Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG

Clonality Polyclonal

Purity Affinity purification


MW 149kDa

Uniprot ID [P19174](#)


Entrez [5335](#)

Dilution Range WB: 1:500-1:2000, IHC-P: 1:10-1:100

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of human placenta



Western blot analysis of extracts of Jurkat cells