



Product Datasheet

ZAP70 (Phospho-Y493) antibody (orb537307)

Description	Rabbit polyclonal antibody to ZAP70 (Phospho-Y493)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic phosphorylated peptide around Y493 of human ZAP70(NP_001070.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70KDa
Uniprot ID	P43403
Dilution Range	WB: 1:2000-1:5000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain...



Western blot analysis of extracts of Jurkat cells...