


## Product Datasheet

### GATC antibody (orb537374)

|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to GATC                                                                                                      |
| <b>Reactivity</b>          | Human                                                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Tested Applications</b> | WB                                                                                                                                      |
| <b>Immunogen</b>           | Recombinant fusion protein of human GATC(NP_789788.1).                                                                                  |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.                                                                                      |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>Isotype</b>             | IgG                                                                                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                                                                              |
| <b>MW</b>                  | 14kDa                                                                                                                                   |
| <b>Uniprot ID</b>          | <b>O43716</b>                                                                                                                           |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                                                                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                         |



Western blot analysis of extracts of MCF...