

Product Datasheet

HUS1 antibody (orb537592)

Description Rabbit polyclonal antibody to HUS1

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, WB

Immunogen Recombinant fusion protein of human HUS1(NP_004498.1).

Preservatives PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 32kDa


Uniprot ID **O60921**

Dilution Range WB: 1:500-1:2000, IHC-P: 1:50-1:200


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of human uterine tissue.



Immunohistochemical staining of human gastric tissue.



Immunohistochemical staining of rat heart tissue.