



Product Datasheet

NPHS1 antibody (orb537603)


Description	Rabbit polyclonal antibody to NPHS1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	A synthetic peptide of human NPHS1
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	130kDa
Uniprot ID	O60500
Entrez	4868
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse kidney...



Immunohistochemical staining of mouse kidney...



Western blot analysis of extracts of various...