


Product Datasheet

DRD4 antibody (orb537647)

Description	Rabbit polyclonal antibody to DRD4
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human DRD4(P21917).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38kDa
Uniprot ID	P21917
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...