


Product Datasheet

PDGFRB (Phospho-Y1021) antibody (orb537894)

Description	Rabbit polyclonal antibody to PDGFRB (Phospho-Y1021)
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y1021 of human Phospho-PDGFRB-Y1021 (NP_002600.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	190kDa
Uniprot ID	P09619
Entrez	5159
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of NIH...