


Product Datasheet

FYN (Phospho-Y416/Lck-Y416) antibody (orb537902)

Description	Rabbit polyclonal antibody to FYN (Phospho-Y416/Lck-Y416)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding Y416 of human SRC Family
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	61kDa
Uniprot ID	P06241
Entrez	53430553932406 , 2 , 023
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of NIH...