



Product Datasheet

ABCB11 antibody (orb538029)


Description	Rabbit polyclonal antibody to ABCB11
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	A synthetic peptide of human ABCB11(NP_003733.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	146kDa
Uniprot ID	O95342
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH/3T3 C...



Immunofluorescence analysis of C6 cells



Western blot analysis of extracts of rat...