


Product Datasheet


PTCH1 antibody (orb538088)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to PTCH1 |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ICC, IF, WB |
| Immunogen | Recombinant fusion protein of human PTCH1(NP_000255.2). |
| Preservatives | PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| Storage | Store at -20°C. Avoid freeze / thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 160kDa |
| Uniprot ID | Q13635 |
| Dilution Range | WB: 1:500-1:2000 |
| Expiration Date | 12 months from date of receipt. |




Immunofluorescence analysis of C6 cells

...



Immunofluorescence analysis of mouse emb...



Immunofluorescence analysis of L-929 cel...