



Product Datasheet

Jak2 (Phospho-Y1007/1008) antibody (orb538262)


Description	Rabbit polyclonal antibody to Jak2 (Phospho-Y1007/1008)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around Y1007 & Y1008 of human JAK2(NP_004963.1?).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	125KDa
Uniprot ID	O60674
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse br...



Immunohistochemical staining of human pl...



Immunohistochemical staining of human lu...