



Product Datasheet

AFG3L2 antibody (orb538364)


Description	Rabbit polyclonal antibody to AFG3L2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Recombinant fusion protein of human AFG3L2(NP_006787.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80-90KDa
Uniprot ID	Q9Y4W6
Dilution Range	WB: 1:200-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ap...



Immunohistochemical staining of human co...



Immunohistochemical staining of human lu...