



Product Datasheet

BCL2L13 antibody (orb538591)

Description	Rabbit polyclonal antibody to BCL2L13
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human BCL2L13(NP_056182.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48kDa
Uniprot ID	Q9BXK5
Dilution Range	WB: 1:500-1:2000, IF: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Western blot analysis of extracts of var...