

Product Datasheet

BRIX1 antibody (orb538629)

Description Rabbit polyclonal antibody to BRIX1

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Recombinant fusion protein of human BRIX1(NP_060791.3).

Preservatives PBS with 0.01% thiomersal,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 41kDa


Uniprot ID **Q8TDN6**

Dilution Range WB: 1:500-1:2000, IHC-P: 1:50-1:100


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of human br...



Immunohistochemical staining of human ap...



Immunohistochemical staining of human lu...