



Product Datasheet

CAMK1 antibody (orb538694)


Description	Rabbit polyclonal antibody to CAMK1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human CaMKI(NP_003647.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41kDa
Uniprot ID	Q14012
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse sp...



Immunohistochemical staining of mouse te...



Immunohistochemical staining of mouse lu...