



Product Datasheet

CHEK2 antibody (orb538871)


Description	Rabbit polyclonal antibody to CHEK2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	Recombinant fusion protein of human Chk2(NP_009125.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65KDa
Uniprot ID	O96017
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100, IF: 1:50-1:200, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of human lu...