



Product Datasheet

CKMT1B antibody (orb538912)


Description	Rabbit polyclonal antibody to CKMT1B
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	Recombinant fusion protein of human CKMT1B(NP_066270.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42kDa
Uniprot ID	P12532
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells



Immunofluorescence analysis of NIH-3T3 cells



Immunohistochemical staining of human colon tissue