

Product Datasheet

DUT antibody (orb539266)

Description Rabbit polyclonal antibody to DUT

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, WB

Immunogen Recombinant fusion protein of human DUT(NP_001020419.1).

Preservatives PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 24kDa


Uniprot ID **P33316**

Dilution Range WB: 1:500-1:1000


Expiration Date 12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells...



Immunohistochemical staining of human appendix...



Immunohistochemical staining of human liver...