

Product Datasheet

EGLN1/EGLN2 antibody (orb539302)

Description Rabbit polyclonal antibody to EGLN1/EGLN2

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen A synthetic peptide of human EGLN1/EGLN2(NP_071334.1/NP_444274.1).

Preservatives PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 44kDa


Uniprot ID [Q96KS0](#), [Q9GZT9](#)

Dilution Range WB: 1:500-1:2000, IHC-P: 1:50-1:200


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of human gastric tissue.



Immunohistochemical staining of human stomach tissue.



Immunohistochemical staining of rat heart tissue.