



Product Datasheet

FLCN antibody (orb539495)


| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to FLCN |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ICC, IF, WB |
| Immunogen | Recombinant fusion protein of human FLCN(NP_659434.2). |
| Preservatives | PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| Storage | Store at -20°C. Avoid freeze / thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 70KDa |
| Uniprot ID | Q8NFG4 |
| Dilution Range | WB: 1:500-1:2000 |
| Expiration Date | 12 months from date of receipt. |



Immunofluorescence analysis of U2OS cell...



Immunofluorescence analysis of NIH/3T3 C...



Immunofluorescence analysis of C6 cells ...