


Product Datasheet

FLVCR2 antibody (orb539503)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to FLVCR2 |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Recombinant fusion protein of human FLVCR2 (NP_060261.2). |
| Preservatives | Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | Affinity purification |
| MW | 57kDa |
| Uniprot ID | Q9UPI3 |
| Entrez | 55640 |
| Dilution Range | WB: 1:500-1:2000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of extracts of var...