



Product Datasheet

GPI antibody (orb539688)


Description	Rabbit polyclonal antibody to GPI
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human glucose-6-phosphate isomerase (NP_000166.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	56kDa
Uniprot ID	P06744
Entrez	2821
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH-3T3 C...



Immunohistochemical staining of rat kidn...



Immunohistochemical staining of rat brai...