



Product Datasheet

HNMT antibody (orb539847)

Description	Rabbit polyclonal antibody to HNMT
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of human HNMT(NP_008826.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33kDa
Uniprot ID	P50135
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse heart tissue.



Western blot analysis of extracts of various tissues.