

## Product Datasheet

### HNRNPD antibody (orb539849)

**Description** Rabbit polyclonal antibody to HNRNPD

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** ICC, IF, IHC, IP, WB

**Immunogen** Recombinant fusion protein of human HNRNPD(NP\_002129.2).

**Preservatives** PBS with 0.05% proclin300,50% glycerol,pH7.3.

**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note** For research use only

**Isotype** IgG


**Clonality** Polyclonal

**MW** 37-45KDa


**Uniprot ID** **Q14103**

**Dilution Range** WB: 1:500-1:2000


**Expiration Date** 12 months from date of receipt.



Immunofluorescence analysis of C6 cells ...



Immunohistochemical staining of mouse brain...



Immunohistochemical staining of human lung...