



Product Datasheet

HPR antibody (orb539871)


| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to HPR |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Recombinant fusion protein of human HPR(NP_066275.3). |
| Preservatives | PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 39kDa |
| Uniprot ID | P00739 |
| Dilution Range | WB: 1:500-1:2000, IHC-P: 1:50-1:100 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemical staining of human pr...



Immunohistochemical staining of mouse lu...



Western blot analysis of extracts of var...