



## Product Datasheet

### IL1B antibody (orb539968)


<b>Description</b>	Rabbit polyclonal antibody to IL1B
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Recombinant fusion protein of human IL1 $\beta$ (NP_000567.1).
<b>Preservatives</b>	PBS with 0.05% proclin300,0.05% rAlbumin,50% glycerol,pH7.3.
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	35KDa/17KDa
<b>Uniprot ID</b>	<b>P01584</b>
<b>Dilution Range</b>	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
<b>Expiration Date</b>	12 months from date of receipt.



Immunofluorescence analysis of mouse lun...



Immunofluorescence analysis of THP-1 cel...



Immunohistochemical staining of rat lung...