


Product Datasheet


ITGB2 antibody (orb540033)

Description	Rabbit polyclonal antibody to ITGB2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human CD18(NP_000202.3).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	95KDa
Uniprot ID	P05107
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of RAW264.7

...



Immunohistochemical staining of Mouse spleen



Immunohistochemical staining of Rat spleen