

## Product Datasheet

### IVNS1ABP antibody (orb540050)

**Description** Rabbit polyclonal antibody to IVNS1ABP

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** WB

**Immunogen** Recombinant fusion protein of human IVNS1ABP(NP\_006460.2).

**Preservatives** PBS with 0.02% sodium azide,50% glycerol,pH7.3.

**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note** For research use only

**Isotype** IgG


**Clonality** Polyclonal

**MW** 72kDa


**Uniprot ID** [Q9Y6Y0](#)

**Dilution Range** WB: 1:500-1:2000, IHC-P: 1:50-1:200


**Expiration Date** 12 months from date of receipt.



Immunohistochemical staining of human to...



Immunohistochemical staining of rat panc...



Immunohistochemical staining of rat brai...