



Product Datasheet

LGMN antibody (orb540221)


Description	Rabbit polyclonal antibody to LGMN
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human Legumain (Legumain (LGMN))(NP_001008530.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36KDa/55KDa
Uniprot ID	Q99538
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of L929 cell...



Immunofluorescence analysis of HeLa cell...



Immunofluorescence analysis of C6 cells

...