

Product Datasheet

MTFP1 antibody (orb540494)

Description Rabbit polyclonal antibody to MTFP1

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Recombinant fusion protein of human MTFP1(NP_057582.2).

Preservatives PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 25kDa


Uniprot ID [Q9UDX5](#)

Dilution Range WB: 1:500-1:2000, IHC-P: 1:50-1:200


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of human br...



Immunohistochemical staining of human CO...



Western blot analysis of extracts of var...