



Product Datasheet

MUL1 antibody (orb540512)


Description	Rabbit polyclonal antibody to MUL1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human MUL1(NP_078820.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45kDa
Uniprot ID	Q969V5
Dilution Range	WB: 1:1000-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of THP-1 cells...



Immunofluorescence analysis of RAW264.7 cells...



Immunohistochemical staining of human liver tissue...