



Product Datasheet

NDUFAB1 antibody (orb540594)


Description	Rabbit polyclonal antibody to NDUFAB1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human NDUFAB1(NP_004994.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	12KDa
Uniprot ID	O14561
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human esophagus tissue.



Immunohistochemical staining of human colon tissue.



Immunohistochemical staining of human colon tissue.