



Product Datasheet

ATF2 (Phospho-T71) antibody (orb540947)

Description	Rabbit polyclonal antibody to ATF2 (Phospho-T71)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around T71 of human ATF2(NP_001871.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70kDa
Uniprot ID	P15336
Dilution Range	WB: 1:500-1:1000, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cel...



Western blot analysis of extracts of var...