

Product Datasheet

PPIG (Phospho-S376) antibody (orb540965)

Description Rabbit polyclonal antibody to PPIG (Phospho-S376)

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen A synthetic phosphorylated peptide around S376 of human PPIG(NP_004783.2).

Preservatives PBS with 0.01% thiomersal,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 110kDa

Uniprot ID [Q13427](#)

Dilution Range WB: 1:500-1:2000

Expiration Date 12 months from date of receipt.



Western blot analysis of extracts of var...