



## Product Datasheet

### PREPL antibody (orb541123)


|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to PREPL                                                                                                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Tested Applications</b> | ICC, IF, WB                                                                                                                             |
| <b>Immunogen</b>           | Recombinant fusion protein of human PREPL(NP_006027.2).                                                                                 |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.                                                                                        |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>Isotype</b>             | IgG                                                                                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                                                                              |
| <b>MW</b>                  | 72kDa                                                                                                                                   |
| <b>Uniprot ID</b>          | <a href="#">Q4J6C6</a>                                                                                                                  |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                                                                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                         |



Immunofluorescence analysis of L929 cell...



Immunofluorescence analysis of C6 cells



Western blot analysis of extracts of var...