



Product Datasheet

PTK2 antibody (orb541193)


Description	Rabbit polyclonal antibody to PTK2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	Recombinant fusion protein of human FAK(NP_722560.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	125Kda
Uniprot ID	Q05397
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200, IP: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of PC12 cell...



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of human st...