

Product Datasheet

RAD51C antibody (orb541259)

Description Rabbit polyclonal antibody to RAD51C

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications WB

Immunogen Recombinant fusion protein of human RAD51C(NP_002867.1).

Preservatives PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 42kDa


Uniprot ID [O43502](#)

Dilution Range WB: 1:500-1:2000, IHC-P: 1:50-1:100


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of mouse st...



Immunohistochemical staining of human st...



Immunohistochemical staining of human pr...