

Product Datasheet

S100A10 antibody (orb541479)

Description Rabbit polyclonal antibody to S100A10

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, WB

Immunogen Recombinant fusion protein of human S100A10(NP_002957.1).

Preservatives PBS with 0.05% proclin300,50% glycerol,pH7.3.

Storage Store at -20°C. Avoid freeze / thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 11KDa


Uniprot ID **P60903**

Dilution Range WB: 1:500-1:2000, IHC-P: 1:50-1:200


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of mouse te...



Immunohistochemical staining of human br...



Immunohistochemical staining of human CO...