



Product Datasheet

ENaC Alpha antibody (orb541516)

Description	Rabbit Polyclonal antibody to ENaC Alpha
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human SCNN1A(NP_001029.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	75KDa
Uniprot ID	P37088
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of mouse kidney



Immunofluorescence analysis of rat lung

...



Immunofluorescence analysis of human lung