



Product Datasheet

SLC23A2 antibody (orb541650)

Description	The absorption of vitamin C into the body and its distribution to organs requires two sodium-depende...
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SLC23A2(NP_005107.4).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70KDa
Uniprot ID	Q9UGH3
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunofluorescence analysis of MCF7 cell...