

## Product Datasheet

### SLC7A2 antibody (orb541712)

**Description** Rabbit polyclonal antibody to SLC7A2

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** IHC, WB

**Immunogen** A synthetic peptide of human SLC7A2(NP\_001008539.3).

**Preservatives** PBS with 0.05% proclin300,50% glycerol,pH7.3.

**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note** For research use only

**Isotype** IgG


**Clonality** Polyclonal

**MW** 72KDa


**Uniprot ID** **P52569**

**Dilution Range** WB: 1:500-1:2000


**Expiration Date** 12 months from date of receipt.



Immunohistochemical staining of mouse te...



Immunohistochemical staining of rat ovar...



Immunohistochemical staining of human pl...