

Product Datasheet

UQCR10 antibody (orb542257)

Description Rabbit polyclonal antibody to UQCR10

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, WB

Immunogen Recombinant fusion protein of human UQCR10(NP_037519.2).

Preservatives PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 10kDa


Uniprot ID [Q9UDW1](#)

Dilution Range WB: 1:500-1:2000, IHC-P: 1:50-1:200


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of rat hear...



Immunohistochemical staining of rat ovar...



Western blot analysis of extracts of var...