


Product Datasheet


USO1 antibody (orb542263)

Description	Rabbit polyclonal antibody to USO1
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 663-962 of human USO1 (NP_003706.2).
Preservatives	PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	108kDa
Uniprot ID	O60763
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of Confocal

...



Immunofluorescence analysis of HeLa cell...



Western blot analysis of extracts of var...