

Product Datasheet

SYNPO2L antibody (orb542730)

Description	Rabbit polyclonal antibody to SYNPO2L
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human SYNPO2L (NP_001107605.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q9H987
Entrez	79933
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.