

## Product Datasheet

### Tafazzin / TAZ antibody (orb542733)

<b>Description</b>	Rabbit polyclonal antibody to Tafazzin / TAZ
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Recombinant fusion protein of human Tafazzin / TAZ (NP_851830.1).
<b>Preservatives</b>	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purification
<b>Uniprot ID</b>	<a href="#">Q16635</a>
<b>Entrez</b>	<a href="#">6901</a>
<b>Dilution Range</b>	WB: 1:1000-1:2000
<b>Expiration Date</b>	12 months from date of receipt.